

Preparation of Rabbit Red Blood Cell Membranes for Transmembrane Pores

Rabbit blood is purchased from www.hemostat.com

1st Wash

- 1a. Take 60 ml raw rabbit blood and divide in 10 ml aliquots in 50 ml centrifuge polypropylene tubes.
- 1b. In each tube, add 40 ml MBSA buffer (150 mM NaCl, 10 mM Na-MOPS, 0.1 % Bovine Serum Albumine (1 mg/ml), titrated with NaOH to reach pH7.4, then filtered with Nalgene Filters (0.2 μ m), and stored at 4 Celsius degrees).
- 1c. This 50 ml solution is centrifugated at 1780 rpm (700 g) for about 10 minutes (in swing bucket).
- 1d. Remove gently the supernatant, and re-suspend gently in MBSA (total volume of 50 ml) buffer again by reversing the 50 ml centrifuge tube upside down.
- 1e. Repeat steps 1c-1d for a total of 4 times.

At the end you expect a 3 ml packed Red blood Cell Membranes

2nd Wash

- 2a Extract packed RBCM. Transfer to new centrifugation tubes.
- 2b. Re-suspend the pellet in lysis buffer up to a final volume of 50 ml (5 mM Na₂HPO₄, 1mM EDTA, pH8.0)
- 2c Centrifugate at high speed centrifuge – 15000 (\geq 20,000 g) for 30 minutes

3rd Wash

- 3a. Collect soft pellet
- 3b. Transfer up to 50 ml centrifuge tubes (equivalent to maxi-prep tubes)
- 3c. Wash the pellet in NaH₂PO₄ until no pink color remain for 10-20 minutes
- 3d. Repeat this procedure about 5 times until no pink color remain

4th Aliquot and freeze - 150 μ l aliquots.

RBLMs were quantified with DC Assay - **Membrane Protein Biorad Protocol**

For DC Assay and Rabbit Red Blood Cell Membranes, please take a look at Steve Cheley et al., Protein Eng. 1997 Dec;10(12):1433-43.